


**LÖWENSTEIN**  
medical



# LUISA

—  
Life support ventilation

With people in mind



# Air for an active life. LUIA.

Inhale, exhale - something most people do without giving it much thought. With LUIA, the new, highly efficient ventilator from Löwenstein Medical, ventilated patients can do the same. LUIA is especially mobile, compact and light-weight. It can be operated in a vertical or horizontal position and flexibly attached to a wheelchair or bed. The ventilator's battery has a run time of up to 18 hours, giving the user a larger radius of action and a good feeling of extra security. It's easy to operate LUIA via the 10-inch touchscreen and intuitive menu navigation. Other features simplify everyday use. LUIA provides the air needed for an active, mobile life and opens up brand new opportunities for users.

- **Battery run time** up to 18 hours
- **Display** 10-inch touchscreen
- **Positioning** vertical or horizontal
- **Connectivity:** LUIA App, prisma therapy software (pTS), monitor connection
- **Monitoring:**  $\text{FiO}_2$  and  $\text{SpO}_2$  measurement
- **Functionality:** ventilation up to 30 ml VT; High-Flow Mode, CPAP and MPV in all circuit systems
- **Safety:** alarm displayed in two languages
- **Uncomplicated:** no adapter, intuitive operation

