


# Popular

## Features & Specifications

- **Construction** : Mild Steel body, Powder Coated for high Corrosion resistance with Anaesthetic appearance.
- **Monitor Tray** : At eye level with Load capacity of 35 kg
- **Vaporiser** : Goldman Halothane Vaporiser and space for 1 No. extra vaporiser.
- **Emergency Press on Block** : Emergency Oxygen cum outlet with pressure release valve facility.
- **Table Top Tray** : Stainless Steel  
Size :427mm x 245mm for Anaesthetic Tools & Ampule Bins for injections & medicines
- **Cylinder Yokes Regulators** : Higher efficiency forged regulators & yokes with S.S. fittings one each for oxygen and nitrous oxide
- **Metal Tubing** : High pressure metal tubing for machine circuit.
- **Mobility** : Four 4" diameter anti-static castor wheels with front brakes.
- **Self Sealing Valve** : Oxygen outlet for Humidifier Bottle. (Single)
- **OFWD** : There will be an audible & Red Light Alarm if the supply of Oxygen fails in Pipeline.
- **Dimension** : H X W X D = 1390mm x 445mm x 490mm audible



Anaesthesia Apparatus  
MADE BY  
Expert Medisystems

Machines are also available in Stainless Steel Trolley

\*For constant upgradation, We take the liberty to change in products & specifications & features without prior notice.